




Above Cuff Vocalization (ACV) for Tracheostomy: Evidence-Based Practice and Clinical Implementation

Carmin Bartow, MS, CCC-SLP



We are passionate about making life easier for people living with a neck stoma, by providing personalized care and innovative solutions.

1

PR0269-06

Disclosures: The following individuals have financial relationship or relationship affiliations to disclose

They are employed by the Educational Division of Atos Medical. There are no other nonfinancial relationships to disclose.

Josie Alston, MA, CCC-SLP

Clinical Education Manager - LARYNGECTOMY & TRACHEOSTOMY

Carmin Bartow, MS, CCC-SLP

Clinical Educator - TRACHEOSTOMY

Megan Quinn, MSN, CPNP-PC

Clinical Educator - TRACHEOSTOMY



2

Disclaimer

The content of this presentation is proprietary and the property of Atos Medical Inc. and is intended solely for the information of the intended recipient. Atos Medical Inc. makes no claim, promise, or guarantee of any kind about the accuracy, completeness, or adequacy of the content of this presentation and expressly disclaims liability for errors and omissions in such content. This presentation may not be shared, copied or distributed without the express written consent of Atos Medical Inc. The content, claims and timelines included in these events, materials and communications are intended for the US market. For information relevant outside of the US, please contact your local Atos Medical office.

Atos

3

Learning Objectives

List the basic steps for using ACV in patients with tracheostomies and mechanical ventilation.

Describe the positive effects of ACV on communication and swallowing.

Identify and troubleshoot potential complications of ACV.

Atos


4

The basics of ACV

TRACOE

ACV is a voicing technique typically used with patients who require cuff inflation during mechanical ventilation

- Tracheostomy tube with a subglottic suction channel is needed
- The suction channel is used to introduce air into the subglottic lumen
- This airflow travels into the upper airway restoring audible voice



Atos

Image: Atos/Tracoe

5

Impact of tracheostomy on communication and swallowing

TRACOE

Atos

6


Mechanical ventilation and the cuff

TRACOE

Purpose of the cuff

To separate the upper and lower airway and prevent air leak during mechanical ventilation

To slow secretions from entering the lower airway




Atos

Image (left): Atos/Tracoe
Image (right): Dräger – Used with permission

7

Let's consider the important functions being bypassed

TRACOE



Atos

Image: Atos/Tracoe

8

Negative impact of tracheostomy on **communication**

Participation, QOL impact

- One of the **most frustrating experiences** reported by patients with tracheostomies and a major factor impacting QOL (Pandian et al. 2014)
- Led to **reduced participation** in treatment (Magnus and Turkington, 2006)
- Restricted the patients' ability to understand information and to participate in care, which led to **universally reported negative emotions**
 - Patients reported feeling of **"frustration", "hopelessness", "confusion", "anger"** related to being unable to communicate (Freeman-Sanderson et al. 2018)

Atos

Pandian et al. (2014); Magnus and Turkington (2006); Freeman-Sanderson et al. (2018)

9

Negative impact of tracheostomy on **communication**

Safety concerns

- Patients with communication difficulties were 3x more likely to experience preventable adverse events than patients without communication problems (Bartlett et al. 2008)

Violation of patient rights

- The Joint Commission, 2010 states that **hospital staff must address the patient's communication needs** before conducting an assessment, providing treatment, obtaining informed consent, discussing end-of-life, or engaging the patient in care discussions

Ineffective pain management

- Difficulty communicating with healthcare providers is a significant barrier to adequate pain relief (Limaye and Katz, 2006)

Atos

Bartlett et al. (2008); Limaye and Katz (2006)

10

Negative impact of tracheostomy on **swallowing**

Does airflow matter?

The mechanoreceptors provide afferent input to the brainstem which modifies the efferent signal for swallowing, key to effective coordinated swallowing

Lack of airflow in the upper airway




Reduced sensory input for swallowing

Atos


11

Implications of a tracheostomy on swallow



Atos

Images: Atos/Tracoe



Raju and Tardi (2021); Yip and Lui (2022); Hess et al. (2014)

12



Reduced sensory input and reduced motor output with tracheostomy

Ponfick et al. (2015) studied patients with critical illness polyneuropathy (CIP) and tracheostomy

- 91% exhibited dysphagia
- 71% exhibited impaired laryngeal sensation

Dysphagia in CIP may result from various contributing factors, including a "learned nonuse" of swallowing muscles

Atos

Ponfick et al. (2015); Robbins et al. (2008)

13



**Therapeutic goal:
Establish upper
airway airflow**

Atos

14

Restoring airflow to the upper airway with tracheostomy **and mechanical ventilation**

TRACOE

- Full or partial cuff deflation for leak speech
- Full cuff deflation and in-line speaking valve
- Above cuff vocalization (ACV)

Atos

15

Leak speech

16

Leak speech

What is it?

- Leaking air around the tracheostomy tube into the upper airway for the purpose of phonation

How to do it?

- SLP / RT teamwork
- Slowly deflate cuff
 - May not need to fully deflate the cuff
- Listen for upper airway sounds / phonation
- Watch for drop in expiratory volumes
- Ventilator adjustments as needed

Atos

17

Leak speech


Suggested readiness parameters

Patient parameters

- Awake, alert, medically stable, adequate airway patency

Vent parameters

- PEEP \leq 10
- FI_{O2} \leq 50%
- PIP \leq 35 - 40



Atos

18

Speaking valve

19

Speaking valve in line with mechanical ventilation

What is it?

- Deflating the cuff
- Placing a PMV in the vent circuit to restore exhalation through the upper airway

How to do it?

- RT / SLP teamwork
- Complete cuff deflation
- If significant drop in expiratory volume, place valve with appropriate adapters
- Ventilator adjustments as needed
- May need continuous monitoring


Atos

20

Cuff deflation and in-line speaking valve

Suggested readiness parameters

- Patient parameters
 - Awake, alert, medically stable, adequate airway patency
- Vent parameters
 - PEEP \leq 10
 - FI_{O2} \leq 50%
 - PIP \leq 35-40



Atos

Image: Atos/Tracoe

21

TRACOE

What if...

The ventilator settings do not meet the readiness criteria for in-line speaking valve?

Your patient and vent settings meet criteria but can only use in-line valve for short periods?

**Above Cuff
Vocalization (ACV)**

Atos

22

Above Cuff Vocalization (ACV)

23


Above Cuff Vocalization

What is it?

- A voicing technique to restore airflow above the inflated cuff and into the subglottic lumen

How to do it?

- Tracheostomy tube with a subglottic suction port is needed
- The suction port is attached to a source of compressed air
- Airflow is introduced into the subglottic lumen
- This airflow travels into the upper airway restoring audible voice



Atos


24

Above Cuff Vocalization

Suggested readiness parameters

Patient parameters

- Awake, alert, medically stable, adequate airway patency



Atos

Image: Atos/Tracoe


25

The basics of ACV

TRACOE

ACV is a voicing technique typically used with patients require cuff inflation during mechanical ventilation

- Tracheostomy tube with a subglottic suction port is needed
- The suction channel is used to introduce air into the subglottic lumen
- This airflow travels into the upper airway restoring audible voice



Atos

Image: Atos/Tracoe

26

TRACOE

Subglottic suction tracheostomy tubes

Atos

27


TRACOE

Tracheostomy tube with subglottic suctioning

Standard tracheostomy tube with the addition of a subglottic suction channel

The additional tubing has small openings seated above the cuff which allows for subglottic suctioning

Can help prevent Ventilator Associated Pneumonia (VAP)



Tracoe Twist Plus Extract

Atos

Image: Atos/Tracoe Papazian LK et al. (2020); Pozuelo-Carrascoso, D. et al. (2020); Muscedere et al. (2011)

28

Tracheostomy tubes with subglottic suction can help prevent ventilator associated pneumonia (VAP)

VAP

Subglottic secretions above the cuff are associated with bacterial colonization of the lower respiratory tract, causing VAP

VAP is a significant cause of morbidity and mortality in critically ill patients who require tracheostomy and mechanical ventilation

Subglottic suctioning

Reduces the volume of secretions

Reduces risk of bacterially contaminated secretions entering the lower respiratory tract

Reduces VAP

For this reason, subglottic suctioning is recommended as one of the preventative measures for VAP


Atos

Terragni et al. (2020); Souza & Santan (2012); Papazian et al. (2020); Pozuelo-Carrascoso et al. (2020); Damas et al. (2015); Muscedere et al. (2011)

29

Subglottic suctioning

The Problem



The Solution




Atos

30

TRACOE

Tracoe Extract tracheostomy tubes feature an innovative design with a flat, subglottic suction channel, and strategically designed and placed suction points. This helps to ensure that secretions are suctioned more rapidly and thoroughly.

Atos



Video: Atos/Tracoe used with permission


31

TRACOE

So, if these trachs can help reduce VAP
AND
 Allow ACV
AND
 Use of in-line speaking valves...
 Why aren't we using them more often??

Great opportunity for RT advocacy

Atos



32

TRACOE

ACV evidence: Communication and swallowing

Atos


33

TRACOE

Above Cuff Vocalization is not new

In the **mid 1960s**, Dr. RML Whitlock described a simple tracheostomy tube attachment to facilitate communication for patients with an inflated cuff.

“The speaking-aid not only makes communication easier but also relieves the patient from the frustration and fear of not being able to make his requirements known”.



Atos

Image: Atos/Tracoe

34

TRACOE


ACV for improved communication

Median time to speech decreased from 13 days to 9 days after implementation of ACV (Gajic et al. 2024)

88% of patients could speak with an audible voice or a whisper (Petosic et al. 2021)

- **Whispered speech is not a fail**

80% of patients achieved an audible voice (McGrath et al. 2018)



Atos
Image: Atos/Tracoe

Gajic et al. (2024); Petosic et al. (2021); McGrath et al. (2018)


35

TRACOE

ACV for improved QOL

ACV improved voice related QOL (V-RQOL) (Petosic et al. 2021)

ACV improved patient-reported QOL (Pandian et al. 2020)



Atos
Image: Atos/Tracoe

Petosic et al. (2021); Pandian et al. (2020)

36

What about delirium?

Communication-vulnerable patients have an increased diagnosis of psychopathology

(JCAHO webinar, *Call to Action: Improving Care to Communication Vulnerable Patients*)



Maybe the inverse is true

Patients who received care through a humanistic lens, including restoration of communication reported enhanced emotional and mental well-being, reduced distress, better understanding of their reality, and decreased delirium. (Batty, 2009)

Atos

Batty (2009)

37


TRACOE

ACV for improved **swallowing, coughing, and secretion management**

Improvements were seen in unstimulated cough, swallow frequency, and aspiration ratings measured during ACV compared to no ACV

Authors postulate improvements may be due to:

- Increase in laryngeal sensitivity and afferent activity
- Physiologic effect of airflow facilitating upward movement of secretions



Atos

Image: Atos/Tracoe

McGrath et al. (2018)

38

ACV for improved **swallowing, coughing, and secretion management**

Kothari et al. (2016) investigated the influence of External Subglottic Airflow (ESAF) on swallowing in patients with severe brain injury

Findings:

Increased swallow frequency during ESAF (ACV)

Reduced volume of secretions after ESAF (ACV)

Atos

Kothari et al. (2016)

39

Tracheostomy Care:
Benefits
of Above Cuff Vocalization



Atos
Learning
Institute

Atos

Video: Atos/Tracoe

40

Candidacy and contraindications

Atos


41

Candidacy considerations

Ideal candidates for ACV when used
for communication:

- Awake and alert
- Have a subglottic suction tracheostomy tube
- Require cuff inflation
 - But not necessarily at all times
- Intact laryngeal function
- No upper airway obstruction
- Have an established stoma (without edema or bleeding)

****Criteria may vary if using ACV to stimulate laryngeal sensation and swallowing***



Atos

Image: Atos/Tracoe

42

TRACOE


Contraindications

Upper airway obstruction

New tracheostomy

- Less than 48 hours post-tracheotomy

Bleeding or edema at stoma



Atos

43

TRACOE

ACV steps and troubleshooting



Atos

44

Multidisciplinary approach

HCPs involved in ACV:

MD
RT
SLP
APN
RN

Staff training is necessary



Atos

45

ACV steps – education

The following ACV steps are not intended to be specific recommendations. They are suggestions for how to approach ACV. You and your teams should create a protocol for ACV at your facility.

Obtain physician order

Explain the procedure to the patient

- Explain what you are doing and what the patient may feel
- Don't guarantee voice



Atos

Image: Atos/Tracoe

46

ACV steps – subglottic suction

Clear the secretions from the subglottic space with a syringe or suction before applying ACV



Atos

Image source: Atos/Tracoe

47

ACV steps – connectors

Connect O2 tubing to the subglottic suction line with or without the thumb port



Without thumb port



With thumb port



Atos

Images: Atos/Tracoe

48

ACV steps – connect O2 tubing to source of compressed air



Atos

Image: Atos/Tracoe

49

ACV steps – introduce airflow

Introduce air slowly through the subglottic suction line into the upper airway starting with 1 LPM

Slowly increase to a typical flow rate of 3–6 LPM depending on patients' needs

Monitor patient's reactions and tolerance

Adjust flow and timing as needed

Flow rates should not exceed 10 LPM

- For many patients this will need to remain lower



Atos

Image: Atos/Tracoe

Petotic et al. (2021)

50

ACV steps – occlude thumb port for voicing and train patient / family (if appropriate and within your facility guidelines)



Atos

Images: Atos/Tracoe

51

TRACOE

Complications and troubleshooting

Discomfort	Endoscopy or FEES for direct visualization
	Head/tube position adjustments
	Reduce flow rate and / or ACV length of time
Dysphonia	Increase flow rate
	Practice / Voice therapy
Air escaping from stoma	Gentle pressure on neck flange
	Adjust vent tubing / vent arms
	Lower flow rate
Emphysema of the neck and face	Immediately discontinue ACV
	Notify MD

Atos

Mills et al. (2022); Calamai et al. (2018)

52

Let's try an experiment...

What is the reported lip-reading accuracy?

- 0-10%
- 10-20%
- 20-40%
- 40-60%

Answer: 10-20%

Based on what we have discussed about effective communication for safety, patient rights, pain management, participation in treatment, QOL, is lip-reading good enough?

Atos

NO

Alteri et al. (2011)

53


Why use ACV?

Atos

54

Benefits of Above Cuff Vocalization

Crucial in delivering patient-centered care
 Does not require cuff deflation
 Does not require alterations to vent settings
 May be able to train patient and family to use ACV
 Improved communication / faster time to voicing
 Promotes laryngopharyngeal sensory stimulation
 Increased swallow frequency
 Reduced secretions
Improvement in QOL



Atos

Image: Atos/Tracoe

Wallace et al. 2022; Kothari, 2017; Petosic et al. 2021; Mills et al. 2022

55





Tracheostomy:

Patient Communication with Above Cuff Vocalization




Video: Atos/Tracoe

56

“Life is so much easier
with a voice” 😊

57

TRACOE

Thank you for
attending!
Any questions?

carmin.bartow@atosmedical.com

References available upon request

Atos

58