

Noninvasive Ventilation

Arthur Jones, EdD, RRT

1

Learning Objectives:

- ^ Explain the rationale for noninvasive ventilation.
- ^ Describe the effects, indications, advantages, disadvantages and complications associated with negative pressure ventilation.
- ^ Describe the operation of specific negative pressure ventilators.
- ^ Describe the modes, effects, complications, indications and contraindications associated with noninvasive positive-pressure ventilation (NPPV).
- ^ Describe the evidence basis for NPPV for selected conditions.

2

Learning Objectives:

- ^ Compare the interfaces used in NPPV with respect to their indications, advantages and limitations.
- ^ Describe the issues pertaining to types of NPPV circuits and humidification systems.
- ^ Compare the ventilator types and modes applied to NPPV with respect to their advantages and disadvantages.
- ^ Describe techniques for ventilator control adjustments for NPPV.
- ^ Discuss clinical issues pertaining to NPPV, including aerosol delivery, heliox, clinical sites and end-of-life care.

3

Negative Pressure Ventilation

4

Definitions

- ▲ Noninvasive ventilation - mechanical ventilation without tracheal tube.
- ▲ Noninvasive positive pressure ventilation (NPPV) - ventilation without tracheal tube and with positive airway pressure
- ▲ Negative pressure ventilation - ventilation with negative pressure applied to thorax.

5


Rationale for NIV

- ▲ Ventilate patients, while avoiding the complications associated with tracheal tubes:
 - ◆ ventilator associated pneumonia
 - ◆ airway trauma
 - ◆ psychological trauma, due to aphonia, restraint

6

NPV physiologic effects

- ▲ Subambient pressure surrounds thorax to inflate lungs



7

Physiologic effects

- ▲ Decreased work of breathing (WOB)
- ▲ Increased distribution of ventilation
- ▲ Subambient pressure inflates lungs
 - ◆ intermittent, with passive deflation
 - ◆ continuous negative expiratory pressure (CNEP) - maintain FRC
- ▲ Increased pulmonary blood flow
- ▲ Increased ventricular filling - increased cardiac output (cuirass)

8

Indications

- ▲ Unable to fit or tolerate mask for NPPV
- ▲ Neuromuscular disease
- ▲ Neurological trauma
- ▲ Intolerance of increased mean airway pressure; e.g., PEEP
- ▲ COPD - chronic state and acute exacerbations

9

Indications

- ▲ Post congenital heart surgery
 - ◆ tetralogy of Fallot correction
 - ◆ tricuspid atresia correction (Fontan)
 - ◆ phrenic nerve injury
- ▲ Neonatal respiratory distress
- ▲ Bronchopulmonary dysplasia

10

Indications

- ▲ Post emphysematous lung resection
- ▲ During microlaryngeal surgery
- ▲ Cystic fibrosis
- ▲ Weaning from PPV
- ▲ Flail chest (CNEP)
- ▲ Meconium aspiration??
- ▲ Bronchopleural fistula??

11

Benefits

- ▲ Avoidance of tracheal tube complications
- ▲ Avoidance of facial trauma from mask
- ▲ Reduced sedation requirements
- ▲ No ventilator - induced lung injury

12

Benefits

- ▲ Patient can talk
- ▲ Patient can cough
- ▲ Improved enteral nutrition - patient may be able to eat

13

Disadvantages

- ▲ Lack of airway protection
- ▲ Large, non-portable equipment (tanks)
- ▲ Decreased patient access (tanks)
- ▲ Cumbersome to apply to some patients (tanks, wraps)

14

Disadvantages

- ▲ Difficult to maintain seal
- ▲ Difficult to monitor volumes
- ▲ Patient intolerance (varies with type)

15

Complications

- ▲ Peripheral venous pooling (tank shock)
- ▲ Gastrointestinal bleeding
- ▲ Dynamic upper airway collapse
- ▲ Irritation at neck seal (tanks)
- ▲ Back pain

16


Contraindications

- ▲ Obstructive sleep apnea
- ▲ Morbid obesity
- ▲ Severe kyphoscoliosis
- ▲ Recent abdominal surgery

17

Negative pressure ventilators

- ▲ Drinker-Shaw iron lung



patient head port

patient bed

18

Negative pressure ventilators

▲ Emerson iron lung



19

Negative pressure ventilators

▲ Porta-Lung™ (tank, only)



20

Negative pressure ventilators

▲ Coppa iron lung and cuirass

- ◆ microprocessor-based
- ◆ not available in the US



21

Negative pressure ventilators

▲ Hayek RTX™



22

Negative pressure ventilators

▲ Pulmo wrap



23

Negative pressure ventilators

▲ Vintage cuirass



24

Negative pressure ventilators

- ▲ Hayek RTX™ cuirass ventilator
 - ◆ Biphasic cuirass ventilation - inspiratory and expiratory pressure
 - ◆ Easy to apply
 - ◆ Modes:
 - f High-frequency chest oscillation
 - f Secretions mode - oscillation , cough
 - f Continuous negative expiratory pressure

Click for video on Hayek ventilator (2.5 min)
https://www.youtube.com/watch?v=t_5EwDP48Pk&t=80s

25

Negative pressure ventilators

- ▲ Hayek RTX™ cuirass ventilator
 - ◆ Modes:
 - f High-frequency chest oscillation
 - f Secretions mode - oscillations, cough
 - f Continuous negative expiratory pressure
 - f Control mode
 - f Patient trigger
 - f Patient synchronized

26

Negative pressure ventilators

- ▲ Hayek RTX™ cuirass ventilator
 - ◆ Specifications
 - f 6-1200 cycles per minute
 - f I:E Ratio: 1:6 - 6:1
 - f Maximum inspiratory pressure: -50 cm H2O

27

Negative pressure ventilators

▲ Hayek RTX™ cuirass ventilator

◆ Specifications

- f6-1200 cycles per minute
- fI:E Ratio: 1:6 - 6:1
- fMaximum inspiratory pressure: -50 cm H₂O
- fMaximum expiratory pressure: +50 cm H₂O
- fPower unit weight: 9 kg
- fFour adult size cuirasses
- fSeven pediatric size cuirasses

28

Ventilator operation

▲ Settings

- ◆ Peak inspiratory pressure - adjusts tidal volume
- ◆ Peak expiratory pressure - active exhalation, cough assistance
- ◆ Continuous negative expiratory pressure
 - f maintains FRC
 - f balances intrinsic PEEP for patient triggering

29

Ventilator operation

▲ Settings

- ◆ Rate
- ◆ I:E ratio
- ◆ Trigger sensitivity
 - f sensed at nares
 - f sensed in cuirass (Hayek)
- ◆ FIO₂ - mask or nasal cannula

30

Monitoring

- ▲ Blood gases
 - ◆ baseline and PRN arterial sampling
 - ◆ pulse oximetry
 - ◆ end-tidal CO₂ monitoring
- ▲ Volumes
 - ◆ spirometry with mask??
 - ◆ respiratory inductive plethysmography??

31

Sites for NPV

- ▲ Intensive care units
- ▲ Intermediate care units
- ▲ Long-term care facilities
- ▲ Homes

32

Patient transport

- ▲ Portable positive pressure ventilation
 - ◆ mouthpiece
 - ◆ mask
- ▲ Iron lung - can be manually operated
- ▲ Battery-powered negative pressure ventilators, like Hayek MRTX™

33

Summary & Review

- ▲ Rationale for NIV - ventilate without intubate
- ▲ Physiologic effects of NPV vs. PPV
 - cardiovascular
- ▲ NPV indications
 - ◆ mask intolerance
 - ◆ need to increase pulmonary perfusion

34

Summary & Review

- ▲ NPV benefits - no ventilator-induced lung injury
- ▲ Complications - tank shock
- ▲ Contraindications - upper airway obstruction
- ▲ NPV enclosure types
 - ◆ tank
 - ◆ cuirass
 - ◆ wrap

35

Summary & Review

- ▲ Ventilator operation
 - ◆ pressure controlled ventilation with supplemental mask or nasal O2
 - ◆ limited monitoring capabilities

36

Noninvasive Positive Pressure Ventilation (NPPV)

37

Attributes

- ▲ Noninvasive positive pressure ventilation - PPV without tracheal tube
- ▲ Important attribute of NPPV - existence of a mask leak - that affects
 - ◆ volume delivered (volume control)
 - ◆ ventilator triggering to inspiration
 - ◆ ventilator cycling to expiration

38

Modes

- ▲ continuous positive airway pressure
- ▲ bilevel positive airway pressure
- ▲ pressure support
- ▲ pressure control
- ▲ volume control
- ▲ proportional assist
- ▲ neurally adjusted ventilatory assist (NAVA)

39

Physiologic effects

- ▲ Decreased WOB
- ▲ Increased dynamic lung compliance
- ▲ Increased tidal volume
- ▲ Increased inspiratory capacity (CPAP and COPD patients)

40

Physiologic effects

- ▲ Improved blood gases
 - ◆ oxygenation increased by end-expiratory pressure
 - ◆ hypercapnea decreased with inspiratory pressure
- ▲ Cardiac output
 - ◆ normal and COPD patients - decreased
 - ◆ some CHF patients - increased

41

Benefits

- ▲ Prevention of ETT complications
- ▲ Reduction in sedation requirements
- ▲ Prevention of tracheotomy
- ▲ Reduction in ICU length-of stay (LOS)

42

Complications

- ▲ Delayed intubation
- ▲ Patient intolerance, anxiety
- ▲ Facial ulcers
- ▲ Ear, sinus pain
- ▲ Increased WOB - patient-ventilator dyssynchrony, due to inappropriate device and/or control settings

43

Complications

- ▲ Pneumothorax
- ▲ Gastric insufflation - high pressures
- ▲ Aspiration
- ▲ Mucus plugging
- ▲ Hemodynamic compromise

44

Contraindications

- ▲ unable to fit or tolerate interface
- ▲ facial trauma or surgery
- ▲ active vomiting
- ▲ acute abdominal process - risk for vomiting, aspiration

45

Contraindications

- ▲ apnea
- ▲ cardiovascular instability
- ▲ excessive and/or viscous secretions
- ▲ recent gastro-oesophageal surgery
- ▲ severely impaired mental status

46

Indications

- ▲ COPD
- ▲ Acute cardiogenic pulmonary edema (ACPE)
- ▲ Blunt thoracic trauma
- ▲ Postoperative respiratory failure
- ▲ Weaning from invasive ventilation

47

Indications

- ▲ COPD
- ▲ Acute cardiogenic pulmonary edema (ACPE)
- ▲ Blunt thoracic trauma
- ▲ Postoperative respiratory failure
- ▲ Weaning from invasive ventilation
- ▲ Miscellaneous conditions
- ▲ Neuromuscular conditions - separate lesson
- ▲ Obstructive sleep apnea - separate lesson

48

Questionable indications

- ▲ ARDS/ALI - may harm by delaying intubation
- ▲ Pneumonia - no evidence of benefit
- ▲ Asthma - no evidence of benefit

49

NPPV & COPD Exacerbations

- ▲ First line treatment for exacerbations
- ▲ Strong evidence for efficacy in hypercapnic failure
- ▲ Effects
 - ◆ decreased WOB
 - ◆ reversal of ventilatory muscle fatigue
 - ◆ decreased PaCO₂
 - ◆ decreased risk for intubation
 - ◆ decreased mortality

50

NPPV & Stable COPD

- ▲ Many COPD patients also have sleep apnea (overlap), with greater risk for hypercapnic failure
- ▲ Effects
 - ◆ Decreased air trapping (TLC)
 - ◆ Increased CO₂ response
 - ◆ Stabilizes heart rhythm by reducing vagal activity

51

NPPV & Stable COPD

^ Situations

- ◆ home - longer survival for adherent patients
- ◆ rehabilitation - may increase exercise tolerance, except for backpack study
- ◆ therapy ceiling for end-stage

52

NPPV for acute cardiogenic pulmonary edema (ACPE)

^ First line treatment - strong evidence

^ Effects - CPAP and bilevel NPPV equally:

- ◆ increased FRC
- ◆ increased lung compliance
- ◆ decreased WOB
- ◆ decreased dyspnea & respiratory rate

53

NPPV for ACPE

^ Effects - CPAP and bilevel NPPV equally:

- ◆ decreased intrapulmonary shunt
- ◆ decreased heart rate
- ◆ increased cardiac output
- ◆ decreased intubation rate
- ◆ may decrease mortality (meta-analysis)

^ Bilevel NPPV may not be more effective than CPAP

54

NPPV & blunt thoracic trauma (flail)

- ^ CPAP and bilevel NPPV studied - fair evidence
 - ◆ case studies
 - ◆ small RCT was stopped due to efficacy of NPPV

55

NPPV & blunt thoracic trauma (flail)

- ^ excluded patients
 - ◆ emergent intubations
 - ◆ injuries to head, face or neck
- ^ effects
 - ◆ decreased rate of pneumonia
 - ◆ decreased mortality

56

NPPV & postop respiratory failure

- ^ CPAP and bilevel NPPV studied
- ^ Postoperative upper abdominal & thoracic surgical patients studied
- ^ Effects:
 - ◆ decreased intubation, reintubation rate
 - ◆ decreased pneumonia, sepsis
 - ◆ decreased mortality
 - ◆ decreased length of hospitalization
- ^ CPAP and bilevel NPPV may be equally effective (need more trials)

57

NPPV & ventilator weaning

▲ Rationales

- ◆ shorten intubation time
 - f decreased sedation
 - f decrease infection
 - f decrease ICU & hospital length-of-stay (LOS)
- ◆ prevent reintubation
- ◆ prevent tracheotomy

58

NPPV & ventilator weaning

▲ Supportive evidence is moderate, when applied to selected patients, who:

- ◆ meet criteria to initiate spontaneous breathing trial
- ◆ meet criteria for extubation:
 - f do not have excessive secretions
 - f have an effective cough
 - f have acceptable mental status
 - f are not a difficult intubation

59

NPPV & ventilator weaning

▲ Supportive evidence is moderate, when applied to selected patients, who:

- ◆ have no impediments for interface
- ◆ tolerate short term spontaneous breathing for mask adjustments, etc.

60

NPPV & post-extubation failure

- ^ Evidence does not support efficacy of NPPV in **treating** post-extubation respiratory failure
- ^ Evidence supports that NPPV may be effective in **preventing** post-extubation respiratory failure where high-risk patients are identified in advance
- ^ NPPV not recommended as a routine intervention for post-extubation situations.

61

NPPV - miscellaneous indications

- ◆ During bronchoscopy to offset increased WOB and hypoxemia
- ◆ Severe bronchiolitis
- ◆ Cystic fibrosis - adults with hypercapneic exacerbations

62

NPPV - miscellaneous indications

- ◆ During bronchoscopy to offset increased WOB and hypoxemia
- ◆ Severe bronchiolitis
- ◆ Cystic fibrosis - adults with hypercapneic exacerbations
- ◆ Immunocompromised patients - prevents ventilator - associated pneumonia
- ◆ Pandemic respiratory infections; e.g., SARS - to prevent infection of caregivers during intubations

63

Summary & Review

- ▲NPPV physiologic effects - may increase cardiac output in CHF
- ▲NPPV benefits - prevent intubation
- ▲NPPV complications - delayed intubation
- ▲NPPV contraindications - mask intolerance
- ▲NPPV indications - cautious application to hypoxemic failure

64

Intermission

65

NPPV Interfaces & Humidification

66

Issues with interfaces

- ▲ Comfort
- ▲ Allowance for patient movement
- ▲ Weight
- ▲ Allergenicity
- ▲ Pressure applied to tissues ==> skin ulceration

67

Issues with interfaces

- ▲ Internal volume
 - ◆ dead space (V_{Drb}) - rebreathed volume
 - ◆ gas compression, decompression volume
- ▲ Leaks - mask seal = 2 cm ==> negligible leaks
- ▲ Multiple sizes available

68

Issues with interfaces

- ▲ Securing system (headgear)
 - ◆ comfort
 - ◆ stability
 - ◆ ease of use
 - ◆ washable for home use
 - ◆ disposable for hospital use

69

Interfaces

- ▲ Mouthpiece
- ▲ Nasal mask
- ▲ Nasal pillows
- ▲ Oronasal mask
- ▲ Total face mask
- ▲ Helmet

70

Mouthpiece interface

- ▲ primarily for daytime use for patients with:
 - ◆ neuromuscular disease
 - ◆ COPD
 - ◆ cystic fibrosis

71

Nasal interfaces

- ▲ Indications
 - ◆ primary interface for obstructive sleep apnea
 - ◆ good starting interface in mild acute respiratory failure, with limitations
 - ◆ postoperative atelectasis (see lung clearance lesson)

72

Nasal interfaces

Advantages:

- ◆ enables speech, eating, coughing
- ◆ less risk for aspiration, gastric distension
- ◆ less claustrophobia

73

Nasal interfaces

Limitations

- ◆ nasal resistance limits effectiveness
- ◆ erroneous monitoring of exhaled TV
- ◆ mouth breathing limits effectiveness and patients with ARF tend to mouth-breathe
- ◆ lesser pressure can be administered

74

Nasal interfaces

masks



75

Nasal interfaces

▲ pillows



76

Oronasal (full face) mask

▲ Most common for bilevel NPPV

▲ Advantages

- ◆ less leakage
- ◆ more stable pressures
- ◆ less patient cooperation

▲ Limitations

- ◆ claustrophobia
- ◆ aspiration

77

Oronasal (full face) mask



78

Total face mask

- ▲ Most effective NPPV interface for acute respiratory failure.
- ▲ Minimal leaks - accommodates greatest pressure
- ▲ Less discomfort
- ▲ Less pressure injury - larger area of contact
- ▲ Does not increase V_{Drb}
- ▲ May increase claustrophobia

79

Total face mask



Respirics Total Face Mask™ (courtesy Respirics)



Respirics Performax™ (courtesy Respirics)

80

Helmet

- ▲ Advantages
 - ◆ researchers compared noninvasive oxygen delivery methods—a helmet versus a face mask—for patients with acute respiratory distress syndrome.
 - ◆ the trial was stopped early because the helmets proved more effective than the face masks for treating this condition.

81

Helmet

^advantages

- ◆ best interface for COVID infection
- ◆ overcomes mask-fit problems
- ◆ more comfortable
- ◆ no facial pressure injury
- ◆ less need for patient cooperation
- ◆ allows speaking, coughing

Click for video on helmet NPPV (2.5 min)

<https://www.youtube.com/watch?v=kuTqecGcwTw>

82

Helmet

^disadvantages

- ◆ not currently FDA-approved
- ◆ may decrease cerebral blood flow in infants
- ◆ impedes patient triggering - decompression volume
- ◆ no capability for volume monitoring
- ◆ humidification may fog the helmet

83

Helmet



84

Helmet



Designer CPAP helmet

Click to hear helmet voice distortion

85

Ventilator circuit

▲ Single-limbed

- ◆ original BiPAP circuit
- ◆ incorporates variable flow leak port
- ◆ requires EPAP > 4 cm H₂O to minimize rebreathing

86

Ventilator circuit

▲ Double-limbed

- ◆ ICU ventilators and recent bilevel ventilators
- ◆ eliminates rebreathing
- ▲ The circuit and interface must be used for the specified ventilator

87

Humidification

- ▲ None - OK for short-term ventilation
humidification ability of mucosa can be overwhelmed
 - ◆ desiccated mucosa releases inflammatory mediators
 - ◆ possible mucus plugging
 - ◆ absence of humidification can impede adherence to therapy

88

Humidification

- ▲ Heat and moisture exchanger
 - ◆ can increase resistance
 - ◆ increase V_{Drb} ==>
 - f hypercapnea
 - f increased minute ventilation
 - f increased WOB

89

Humidification

- ▲ Heat and moisture exchanger
 - ◆ increase V_{Drb}
 - ◆ can increase resistance
- ▲ Heated humidification
 - ◆ no effect on V_{Drb}
 - ◆ no effect on resistance
- ▲ Ambient temperature passive humidification increases comfort for some patients

90

Summary & Review

- ▲ Issues with NPPV interfaces - comfort, pressure injury, leaks
- ▲ Specific interfaces:
 - ◆ mouthpiece
 - ◆ nasal mask, pillow
 - ◆ oronasal mask
 - ◆ total face mask
 - ◆ helmet - not FDA-approved

91

Summary & Review

- ▲ NPPV circuits - single, vs. double-limbed
- ▲ Humidification - HMEs and V_{Drb}

92

NPPV Devices & Controls

93

Institutional bilevel ventilators

▲ Desirable features

- ◆ built-in blender
- ◆ leak compensation
- ◆ trigger compensation
- ◆ backup rate
- ◆ rise time adjustment
- ◆ graphic display
- ◆ alarms
- ◆ battery power supply

94

Institutional bilevel ventilators

▲ Resironics

- ◆ BiPAP S/T™ - no blender
- ◆ Focus™ - no blender
- ◆ Vision™
- ◆ V60

95

Institutional bilevel ventilators

▲ Resironics Vision™

- ◆ acute care bilevel ventilation
- ◆ blender - adjustable FIO2
- ◆ Auto-Track™ trigger
- ◆ rise time adjustable
- ◆ graphic display
- ◆ alarms
- ◆ no battery backup

96

Institutional bilevel ventilators

▲ Respironics V60™

- ◆ acute care bilevel ventilation
- ◆ average volume-assured pressure support
- ◆ auto-adaptive leak compensation
- ◆ inspiratory triggering
- ◆ expiratory cycling
- ◆ adjustable ramp time
- ◆ built-in blender
- ◆ battery backup- 6 hours

97

ICU ventilators

▲ With NIV mode - provide leak compensation

▲ Without NIV mode - no leak compensation

- ◆ impedes patient triggering
- ◆ impedes cycling to expiration in PSV mode
- ◆ can be used; but, could lead to problems (injury, litigation)

98

Ventilators

▲ ICU ventilators with noninvasive modes (examples)

- ◆ Maquet Servo™
- ◆ Drager Evita XL™
- ◆ Newport e360™ and HT50™ (transport)
- ◆ Viasys Vela™ and Avea™
- ◆ Hamilton Raphael™ & C-2™
- ◆ GE Engstrom Carestation™
- ◆ Puritan-Bennett 840™

99

Ventilators

▲ Ventilator Intervention Technology Accessible Locally (VITAL)

- ◆ developed by NASA for pandemic
- ◆ tested in Jet Propulsion Laboratory (JPL)
- ◆ FDA emergency authorization April 30, 2020
- ◆ uses some parts from available ventilators
- ◆ eight U.S. manufacturers licensed for production.

100

Ventilators

▲ Ventilator Intervention Technology Accessible Locally (VITAL)

- ◆ assist-control
- ◆ volume targeted
- ◆ PEEP 5-35 cm H₂O
- ◆ frequency 4-40 BPM
- ◆ TV 150-800 mL

101

Ventilators

▲ Ventilator Intervention Technology Accessible Locally (VITAL)



102

ICU ventilators


▲ advantages

◆ graphics to show:

- ftriggering
- f Cycling
- f intrinsic PEEP (PEEPi)
- f work of breathing (WOB)

103

Spontaneous WOB



104

ICU ventilators

▲ advantages

- ◆ graphics
- ◆ exhaled volume monitoring
- ◆ expiratory cycle adjustment
- ◆ rise time adjustment
- ◆ invasive ventilation capability - easier to switch over

▲ disadvantage - expense

105

Modes applied to NIV

▲ volume control

- ◆ identical success rate, compared to pressure control
- ◆ stable volumes in face of changing lung mechanics
- ◆ higher peak airway pressures that caused flatulence in two patients that sounded like a mask leak
- ◆ common in home NIV (Europe)

106

Modes applied to NIV

▲ pressure control ventilation (PCV)

- ◆ stable volume delivery in face of leaks
- ◆ flow variable with patient demands
- ◆ trial results - PCV may be more effective than PSV for COPD exacerbations

107

Modes applied to NIV

▲ pressure support

- ◆ most common mode for NIV
- ◆ some ventilators do not include backup rate - apnea adjustment is important

108

Modes applied to NIV

- ▲ pressure support
 - ◆ appropriate expiratory cycle adjustment is important
 - ◆ appropriate rise time adjustment is important
 - ◆ observe inspiratory time and I:E ratio - may require inspiratory time limit adjustment

109

Modes applied to NIV

- ▲ proportional assist ventilation (PAV)
 - ◆ delivers flow proportional to patient's inspiratory effort
 - ◆ terminates flow in response to cessation of inspiratory effort

110

Modes applied to NIV

- ▲ proportional assist ventilation (PAV)
 - ◆ clinical trials conclude that it is better tolerated than PSV; but, PSV was delivered with PB 7200ae?
 - ◆ Respirationics Vision - PAV mode not available in U.S.
 - ◆ Puritan-Bennett 840 - company does not support use of PAV with NIV

111

Ventilation modes and NIV

- ▲ Pressure support with volume guarantee (VSV)??
 - ◆ no NIV with VSV trials located
 - ◆ problem with VSV - patient distress and hyperpnea causes ventilator to decrease support

112

Ventilation modes and NIV

- ▲ Neurally adjusted ventilatory assist (Maquet NAVA™)
 - ◆ flow delivery in response to diaphragmatic electrical activity
 - ◆ eliminates leaks as a triggering and cycling factor
 - ◆ evidence favors NAVA over PSV

113

Ventilation settings

- ▲ PEEP 3 < > 10 (ideally)
- ▲ PSV, IPAP - exhaled tidal volume > 5 ml/kg IBW
- ▲ Reasonable starting pressures - PEEP = 5; PSV = 10
- ▲ FIO₂ for desired SPO₂

114

Ventilation settings

▲ Expiratory flow cycling adjustment:

- ◆ patient comfort
- ◆ ventilator graphics
- ◆ inspiratory time, I:E ratio
- ◆ to eliminate PEEP_i


115

Expiratory flow-cycling (PSV)

Patient expiratory effort

Late termination

Inability to trigger



Courtesy Newport Medical

116

Ventilation settings

▲ Rise time adjustment:


- ◆ patient comfort
- ◆ ventilator graphics

117

Inspiratory flow/rise- pressure wave

Linear or bowed upward rise in pressure after trigger on the pressure wave

Slow rise in pressure, concave shape of the pressure wave



Courtesy Newport Medical

118

Summary & Review

- ▲ Types of ventilators used for NPPV
- ▲ Specific bilevel and ICU ventilators
- ▲ ICU ventilator advantages
- ▲ NPPV modes
- ▲ Ventilation settings

119

NPPV Clinical Issues

120

Clinical indicators of successful NIV

- ▲ Favorable response in first 2 hours of NPPV
- ▲ APACHE II score < 29
- ▲ pH > 7.30
- ▲ Glasgow coma score \geq 15

121

Clinical indicators of successful NIV

- ▲ Favorable response in first 2 hours of NPPV
- ▲ APACHE II score < 29
- ▲ pH > 7.30
- ▲ Glasgow coma score \geq 15
- ▲ Absence of pneumonia or ARDS
- ▲ Minimal interface air leaks
- ▲ Patient-ventilator synchrony

122

Aerosol delivery

- ▲ Remove NIV interface and administer, if safe and feasible
- ▲ Aerosol medications are effective with NIV via all types of devices
- ▲ Aerosols are effectively delivered via nasal interface
- ▲ Increased dosage may be necessary, because of leaks, nasal deposition

123

Aerosol delivery

- ▲ Aerosol delivery devices
 - ◆ pneumatic nebulizer
 - ∫ place proximal to patient
 - ∫ place between patient and leak valve
 - ◆ vibrating mesh (Aeroneb Pro™) - does not add flow to circuit
 - ◆ MDI & spacer - coordinate with inspiration

124

Aerosol delivery

- ▲ Total face masks allow medications to enter eyes - especially problematic with ipratropium
- ▲ More research is needed to generate specific recommendations
- ▲ Regardless of device, it's a good idea to place aerosol generator between the HME and the patient

125

Heliox therapy

- ▲ Meta-analyses on routine use of heliox for asthma and COPD do not support routine use
- ▲ Heliox decreases PEEPi and WOB in ventilated patients with COPD exacerbations

126

Heliox therapy

▲ Problems:

- ◆ heliox can not be used in hypoxemic patients
- ◆ heliox causes errors in TV and FIO2 measurement and delivery
- ◆ only FDA-approved ventilator is the Viasys Avea
- ◆ heliox is expensive
- ◆ helium causes vocal distortion

127

Sites for NPPV administration

▲ Factors

- ◆ acuity of patient
- ◆ expertise of personnel
- ◆ availability of physical resources

▲ Monitoring capabilities

- ◆ personnel (skilled)
- ◆ electronic

128

Sites for NPPV administration

▲ Pre-hospital - EMS

- ◆ for some conditions, the sooner, the better for NPPV
- ◆ avoids emergency intubations
- ◆ especially applicable to ACCPE & COPD

▲ Emergency room

- ◆ early initiation of NPPV
- ◆ advanced resources, including RTs

129

Sites for NPPV administration

- ▲ Intensive care
 - ◆ best site for sickest patients
 - ◆ intensive monitoring
 - ◆ extensive physical resources
 - ◆ personnel resources
 - ↳ respiratory therapists
 - ↳ critical care nurses

130

Sites for NPPV administration

- ▲ Intermediate care (step-down)
 - ◆ usually telemetric monitoring
 - ◆ personnel resources varies:
 - ↳ respiratory therapists
 - ↳ patient:nurse ratio
 - ◆ stable patients

131

Sites for NPPV administration

- ▲ General ward
 - ◆ telemetry - maybe
 - ◆ more patients per nurses, who may be unfamiliar with NPPV
 - ◆ respiratory therapy coverage varies
 - ◆ intermittent NPPV, as for:
 - ↳ sleep apnea
 - ↳ stable COPD
 - ↳ stable neuromuscular disease

132

Sites for NPPV administration

- ▲ Long-term care facilities
 - ◆ chronic care
 - f COPD
 - f failure to wean from ventilation
 - ◆ monitoring varies by units
 - ◆ usually have skilled respiratory therapy staff

133

Sites for NPPV administration

- ▲ Home
 - ◆ chronic conditions
 - ◆ end-of-life care
 - ◆ requires education of patient & caregivers

134

NIPPV and End-of-Life Care

- ▲ Patient choices
 - ◆ do not intubate (DNI)
 - ◆ comfort measures only (CMO)
- ▲ Informed consent of patient and/or family is needed - NIV is life support
- ▲ Common conditions
 - ◆ COPD
 - ◆ cancer
 - ◆ neuromuscular diseases
 - ◆ chronic heart failure

135

NIPPV and End-of-Life Care

▲ Goals of NIV for terminal patients

◆ delay death

- f to go home
- f to settle personal issues
- f to see a person

◆ provide comfort - for whom?

- f decrease dyspnea
- f comfort is not provided when a patient is resisting the treatment

136

NIPPV and End-of-Life Care

▲ Ethical controversy exists over whether NIV ought to be used at end-of-life.

▲ The decision should rest with the patient (author's opinion)

137

Summary & Review

- ▲ Indicators for successful NPPV
- ▲ Aerosol delivery - does work with NPPV
- ▲ Heliox - may decrease WOB for COPD
- ▲ Sites for NPPV delivery
- ▲ NPPV for end-of-life care

138

END

139

References

- ▲ Sweet DD, Naismith A, Keenan SP, Sinuff T, Dodek PM. Missed opportunities for noninvasive positive pressure ventilation: a utilization review. *J Crit Care.* 2008 Mar;23(1):111-7. Epub 2007 Jul 5.
- ▲ Budweiser S, Jörres RA, Pfeifer M. Noninvasive home ventilation for chronic obstructive pulmonary disease: indications, utility and outcome. *Curr Opin Pulm Med.* 2008 Mar;14(2):128-34.
- ▲ Corrado A, Gorini M. Negative pressure ventilation: Is there still a role? *Respir J* 2002;20:187-197.
- ▲ Corrado, A, et al. Iron lung versus conventional mechanical ventilation in acute exacerbation of COPD. *Eur Respir J.* 2004 Mar ;23 (3):419-24.
- ▲ Corrado, A, et al. Negative pressure ventilation versus conventional mechanical ventilation in the treatment of acute respiratory failure in COPD patients. *Eur Respir J* 1998 12: 519-525

140

References

- ▲ Shekerdeman LS., Bush A, Shore DF, Lincoln C, Redington AN. Cardiorespiratory responses to negative pressure ventilation after tetralogy of Fallot repair: a hemodynamic tool for patients with a low-output state *J Am Coll Cardiol* 1999 33: 549-555
- ▲ Wedzicha, J.A., Muir, J-F. Noninvasive ventilation in chronic obstructive pulmonary disease, bronchiectasis and cystic fibrosis. *Eur Respir J* 2002 20: 777-784.
- ▲ Kössler W, Lahrmann H, Brath H, Wei T, Frank W, Wild M, Zwick H, Wanke T. Feedback-Controlled Negative Pressure Ventilation in Patients with Stable Severe Hypercapnic Chronic Obstructive Pulmonary Disease. *Respiration* 2000;67:362-366.

141

References

- ▲ Smith, IE, King, MA, Shneerson, JM. Choosing a negative pressure ventilation pump: are there any important differences? *Eur Respir J* 1995 8: 1792-1795.
- ▲ Zibrak JD, Hill NS, Federman EC, Kwa SL, O'Donnell C. Evaluation of intermittent long-term negative-pressure ventilation in patients with severe chronic obstructive pulmonary disease. *Am Rev Respir Dis*. 1988 Dec;138(6):1515-8.
- ▲ Corrado A, Gorini G, Vilella G, De Paolua E. Negative pressure ventilation in the treatment of acute respiratory failure: an old noninvasive technique reconsidered. *Eur Respir J* 1996;9:1531-1544.

142

References

- ▲ Kallet RH. Noninvasive ventilation in acute care: controversies and emerging concepts. *Respir Care*. 2009 Feb;54(2):259-63.
- ▲ Gay PC. Complications of noninvasive ventilation in acute care. *Respir Care*. 2009 Feb;54(2):246-58.
- ▲ Kacmarek RM. Should noninvasive ventilation be used with the do-not-intubate patient? *Respir Care*. 2009 Feb;54(2):223-31.
- ▲ Epstein SK. Noninvasive ventilation to shorten the duration of mechanical ventilation. *Respir Care*. 2009 Feb;54(2):198-211.
- ▲ Mehta S, Al-Hashim AH, Keenan SP. Noninvasive ventilation in patients with acute cardiogenic pulmonary edema. *Respir Care*. 2009 Feb;54(2):186-97.

143

References

- ▲ Keenan SP, Mehta S. Noninvasive ventilation for patients presenting with acute respiratory failure: the randomized controlled trials. *Respir Care*. 2009 Jan;54(1):116-26.
- ▲ Kallet RH, Diaz JV. The physiologic effects of noninvasive ventilation. *Respir Care*. 2009 Jan;54(1):102-15.
- ▲ Nava S, Navalesi P, Gregoret C. Interfaces and humidification for noninvasive mechanical ventilation. *Respir Care*. 2009 Jan;54(1):71-84.
- ▲ Davies JD, Gentile MA. What does it take to have a successful noninvasive ventilation program? *Respir Care*. 2009 Jan;54(1):53-61.
- ▲ Scala R, Naldi M. Ventilators for noninvasive ventilation to treat acute respiratory failure. *Respir Care* 2008;53:1054-1080.

144

References

- ▲ Xirouchaki N, Kondoudaki E, Anastasaki M, Alexopoulou C, Koumiotaki S, Georgopoulos D. Noninvasive bilevel positive pressure ventilation in patients with blunt thoracic trauma. *Respiration*. 2005 Sep-Oct;72(5):517-22. Epub 2005 Jul 1.
- ▲ deBoisblanc MW, Goldman RK, Mayberry JC, Brand DM, Pangburn PD, Soifer BE, Mullins RJ. Weaning injured patients with prolonged pulmonary failure from mechanical ventilation in a non-intensive care unit setting. *J Trauma*. 2000 Aug;49(2):224-30;
- ▲ Bach JR, Alba AS. Noninvasive options for ventilatory support of the traumatic high level quadriplegic patient. *Chest*. 1990 Sep;98(3):613-9.
- ▲ Linton DM, Sviri S. Conservative management of flail chest after cardiopulmonary resuscitation by continuous negative extrathoracic pressure. *Crit Care Resusc*. 2006 Dec;8(4):339-40.

145

References

- ▲ Al-Ansari MA. Successful use of non-invasive positive pressure ventilation in a complicated flail chest. *Saudi Med J*. 2006 Aug;27(8):1244-7.
- ▲ Gunduz M, Unlugenc H, Ozalevli M, Inanoglu K, Akman H. A comparative study of continuous positive airway pressure (CPAP) and intermittent positive pressure ventilation (IPPV) in patients with flail chest. *Emerg Med J*. 2005 May;22(5):325-9.
- ▲ Van de Louw A, Brocas E, Boiteau R, Perrin-Gachadoat D, Tenailon A. Esophageal perforation associated with noninvasive ventilation: a case report. *Chest*. 2002 Nov;122(5):1857-8.
- ▲ Girault C, Briel A, Benichou J, Hellot MF, Dachraoui F, Tamion F, Bonmarchand G. Interface strategy during noninvasive positive pressure ventilation for hypercapnic acute respiratory failure. *Crit Care Med*. 2009 Jan;37(1):124-31.

146

References

- ▲ De Santo LS, Bancone C, Santarpino G, Romano G, Della Corte A, Vicchio M, De Pietro A, Galdieri N, Cotrufo M. Noninvasive positive-pressure ventilation for extubation failure after cardiac surgery: Pilot safety evaluation. *J Thorac Cardiovasc Surg*. 2009 Feb;137(2):342-6.
- ▲ Nickol AH, Hart N, Hopkinson NS, Hamnegård CH, Moxham J, Simonds A, Polkey MI. Mechanisms of improvement of respiratory failure in patients with COPD treated with NIV. *Int J Chron Obstruct Pulmon Dis*. 2008;3(3):453-62.
- ▲ Tomii K, Seo R, Tachikawa R, Harada Y, Murase K, Kaji R, Takeshima Y, Hayashi M, Nishimura T, Ishihara K. Impact of noninvasive ventilation (NIV) trial for various types of acute respiratory failure in the emergency ICU department; decreased mortality and use of the ICU. *Respir Med*. 2009 Jan;103(1):67-73.

147

References

- ▲ Lin MS, Guo HR, Huang MH, Chen CR, Wu CL. Predictors of successful noninvasive ventilation treatment for patients suffering acute respiratory failure. *J Chin Med Assoc.* 2008 Aug;71(8):392-8.
- ▲ Creagh-Brown B, Shee C. Noninvasive ventilation as ceiling of therapy in end-stage chronic obstructive pulmonary disease. *Chron Respir Dis.* 2008;5(3):143-8.
- ▲ Nava S, Cirio S, Fanfulla F, Carlucci A, Navarra A, Negri A, Ceriana P. Comparison of two humidification systems for long-term noninvasive mechanical ventilation. *Eur Respir J.* 2008 Aug;32(2):460-4.
- ▲ Scala R, Naldi M. Ventilators for noninvasive ventilation to treat acute respiratory failure. *Respir Care.* 2008 Aug;53(8):1054-80.

148

References

- ▲ Vital FM, Saconato H, Ladeira MT, Sen A, Hawkes CA, Soares B, Burns KE, Atallah AN. Non-invasive positive pressure ventilation (CPAP or bilevel NPPV) for cardiogenic pulmonary edema. *Cochrane Database Syst Rev.* 2008 Jul 16;(3):CD005351.
- ▲ Khayat RN, Abraham WT, Patt B, Roy M, Hua K, Jarjoura D. Cardiac effects of continuous and bilevel positive airway pressure for patients with heart failure and obstructive sleep apnea: a pilot study. *Chest.* 2008 Dec;134(6):1162-8. Epub 2008 Jul 18.
- ▲ Gray A, Goodacre S, Newby DE, Masson M, Sampson F, Nicholl J. Noninvasive ventilation in acute cardiogenic pulmonary edema. *N Engl J Med.* 2008 Jul 10;359(2):142-51.

149

References

- ▲ Antonelli M, Bello G. Noninvasive mechanical ventilation during the weaning process: facilitative, curative, or preventive? *Crit Care.* 2008;12(2):136. Epub 2008 Apr 21.
- ▲ Trevisan CE, Vieira SR; Research Group in Mechanical Ventilation Weaning. Noninvasive mechanical ventilation may be useful in treating patients who fail weaning from invasive mechanical ventilation: a randomized clinical trial. *Crit Care.* 2008;12(2):R51. Epub 2008 Apr 17.
- ▲ Bellani G, Patroniti N, Greco M, Foti G, Pesenti A. The use of helmets to deliver non-invasive continuous positive airway pressure in hypoxemic acute respiratory failure. *Minerva Anesthesiol.* 2008 Nov;74(11):651-6. Epub 2008 Jul 18.

150

References

- ^ Dreher M, Doncheva E, Schwoerer A, Walterspacher S, Sonntag F, Kabitz HJ, Windisch W. Preserving Oxygenation during Walking in Severe Chronic Obstructive Pulmonary Disease: Noninvasive Ventilation versus Oxygen Therapy. *Respiration*. 2008 Dec 18.
- ^ Budweiser S, Hitzl AP, Jörres RA, Heinemann F, Arzt M, Schroll S, Pfeifer M. Impact of noninvasive home ventilation on long-term survival in chronic hypercapnic COPD: a prospective observational study. *Int J Clin Pract*. 2007 Sep;61(9):1516-22.
- ^ Victor J, Moran JL, Phillips-Hughes J, Graham P, Bersten AD. Effect of noninvasive positive pressure ventilation on mortality in patients with acute cardiogenic pulmonary edema: a meta-analysis. *Lancet* 2006;367:1155-63.

151

References

- ^ Burns KE, Adhikari NK, Keenan SP, Meade M. Use of non-invasive ventilation to wean critically ill adults off invasive ventilation: meta-analysis and systematic review. *BMJ* 2009; 338.b1574.
<http://www.bmj.com/content/338/bmj.b1574.full.pdf>
- ^ Holanda MA, Reis RC, Winkeler GF, Fortaleza SC, Lima JW, Pereira ED. Influence of total face, facial and nasal masks on short-term adverse effects during noninvasive ventilation. *J Bras Pneumol*. 2009 Feb;35(2):164-73.
- ^ Rana S, Jenad H, Gay PC, Buck CF, Hubmayr RD, Gajic O. Failure of non-invasive ventilation in patients with acute lung injury: observational cohort study. *Crit Care*. 2006;10(3):R79. Epub 2006 May 12.

152