

# **Pulmonary Function Testing Part I**

Arthur Jones, EdD, RRT

# **Introduction to Pulmonary Function Testing**

## **Learning Objectives:**

- ^ Describe the devices and methods for direct/ indirect spirometry and body plethysmography, for lung volume and compliance measurement.
- ^ Interpret lung volumes and compliance
- ^ Describe the devices and methods for measuring inspiratory and expiratory pressures.
- ^ Interpret inspiratory and expiratory pressures.

## **Purposes of PFTs:**

- ^ Suggest a diagnosis by identifying the pattern of disease; e.g.,
  - ◆ restrictive
  - ◆ obstructive
  - ◆ neuromuscular
  - ◆ combined defects
- ^ Assess the severity and progression of disease

## **Learning Objectives:**

- ^ Describe the physiologic basis and indications for spirometric flow measurements.
- ^ Describe the devices and methods for spirometric flow measurements.
- ^ Interpret results from spirometric flow measurements.
- ^ Describe body plethysmographic measurement of airway resistance.

## **Purposes of PFTs:**

- ^ Evaluate therapeutic & rehabilitation regimens
- ^ Evaluate surgical risk
- ^ Determine disability
- ^ Estimate prognosis

### Types of PFTs:

- ^Lung volumes
  - ◆direct spirometry
  - ◆dilution methods
  - ◆plethysmography
- ^Lung mechanics
  - ◆airway flow
  - ◆airway resistance
  - ◆lung compliance
- ^Diffusing capacity


### Lung compartments

- ^Lung volumes
  - ◆Tidal volume (TV)
  - ◆Inspiratory reserve volume (IRV)
  - ◆Expiratory reserve volume (ERV)
  - ◆Residual volume (RV)

### Types of PFTs:

- ^Bronchodilator response
- ^Bronchoprovocation
- ^Exercise testing
- ^Metabolic testing

### Lung compartments




AP Jones

## Lung Volume Measurement

### Lung compartments

- ^Lung capacities two or more volumes
  - ◆Inspiratory capacity (IC) = TV + IRV
  - ◆Functional residual capacity (FRC) = ERV + RV
  - ◆Vital capacity (VC) = IRV + TV + ERV
  - ◆Total lung capacity (TLC) = TV + IRV + ERV + RV

## Lung compartments



## Direct spirometry devices

- ▲ volume displacement spirometers:
  - ◆ water seal; e.g., Collins, Tissot




## Direct spirometry

- ▲ Measures all volumes and capacities that do NOT include the RV; i.e., FRC & TLC

## Direct spirometry


- ▲ volume displacement spirometers
  - ◆ dry rolling seal; e.g., Spirotech S780, nSpire HD CPL



AP  
Jones

## Direct spirometry devices

- ▲ volume displacement spirometers:
  - ◆ water seal; e.g., Collins, Tissot



## Direct spirometry


- ▲ volume displacement spirometers
  - ◆ dry rolling seal; e.g., Spirotech S780



## Direct spirometry

### ^ volume displacement spirometers

#### ◆ bellows



## Direct spirometry

### ^ volume displacement spirometers

### ^ mechanical respirometers

### ^ electronic spirometers - flow sensors

#### ◆ pneumotachometers

#### ◆ heated wire sensors

#### ◆ pitot tube sensors

#### ◆ ultrasonic flow sensors

## Direct spirometry

### ^ volume displacement spirometers

#### ◆ bellows - Vitalograph Gold Standard™



## Indirect spirometry

### ^ Gas dilution methods

#### ◆ measure capacities that contain RV

#### ◆ Principle:

$$C_1V_1 = C_2V_2 \implies$$

#### ◆ $V_1(\text{FRC}) = C_2V_2/C_1$

## Direct spirometry

### ^ Volume displacement spirometers

### ^ Turbines - flow sensors; e.g., Wright respirometer



## Gas dilution methods

### ^ Methods

#### ◆ Helium dilution

#### ◆ Nitrogen washout

## Nitrogen washout

- △ Open circuit method
- △ Patient breathes 100% O<sub>2</sub>, starting at FRC
- △ Exhaled gas is collected in Tissot spirometer or Douglas bag
- △ N<sub>2</sub> concentration measured
- △ When minimal N<sub>2</sub>% (1%) is reached, volume & %N<sub>2</sub> in Tissot are measured

## Nitrogen washout


- △ Calculation includes correction factors

$$FRC = \frac{(V_{exh})(\%N_{2,exh})}{\%N_{2,alv} (.75)}$$

## Nitrogen washout

- △ Douglas collection bag

Typical Douglas Bag Method Components



## Helium dilution method

- △ Closed circuit method
- △ Patient rebreathes HeO<sub>2</sub> mixture with known volume of He at FRC until equilibration between lung and spirometer
- △ O<sub>2</sub> is added during the test
- △ CO<sub>2</sub> is absorbed during the test

## Nitrogen washout

- △ Tissot water seal spirometer



## Helium dilution method

- △ Total volume for circuit and spirometer calculated (with correction factors)

$$FRC = \frac{(\%He_{initial} - \%He_{final}) \times \text{system volume}}{\%He_{final}}$$

### Helium dilution method

#### ^ Sources of error

- ◆ patient intolerance - long procedure
- ◆ leaks in system
- ◆ leaks in patient; e.g., eardrums
- ◆ slow, or non-communicating lung units
  - slow units increase equilibration time
  - closed units cause underestimation of FRC

### Whole body plethysmography

#### ^ Advantages

- ◆ measures all thoracic gas
- ◆ rapid measurement
- ◆ better cooperation
- ◆ also measures  $R_{AW}$

### Whole body plethysmography

^ Measures thoracic gas volume (TGV,  $V_{TG}$ ) - includes all gas in the thorax

^ Trapped gas can be estimated:

$$\text{Volume}_{\text{trapped}} = \text{TGV}_{\text{pleth}} - \text{FRC}_{\text{He dilution}}$$

^ Principle - Boyle's law

- ◆  $P_1V_1 = P_2V_2$

### Whole body plethysmography

#### ^ Disadvantages

- ◆ claustrophobia
- ◆ equipment expense
- ◆ equipment space requirement
- ◆ extensive operator training
- ◆ may overestimate FRC

### Whole body plethysmography

- ^ Patient sits in sealed chamber
- ^ Pants at FRC against closed shutter
- ^ Mouth pressure change is measured
- ^ Chamber volume change is measured
- ^ Oscilloscope displays relationship of changes in pressure & volume (tangent)

- ◆  $\text{TGV } (V_{TG}) = (dV/dP_A) (P_A - dP_A)$

### Alternate methods for lung volumes

#### ^ Imaging techniques

- ◆ standard radiographs
- ◆ computed tomography
- ◆ magnetic resonance imaging

^ Oxygen washout - during mechanical ventilation

### Abnormal volumes - implications

- △ Reduced volumes with decreased lung compliance ( $C_L$ ) ==> restriction
  - ◆ intrinsic pulmonary restrictive disease
    - ▶ fibrotic processes; e.g., silicosis, idiopathic pulmonary fibrosis
    - ▶ pulmonary edema
    - ▶ space occupying lesions; e.g., tumors
    - ▶ surfactant deficiency

### Abnormal volumes - implications

- △ Reduced volumes
  - ◆ additional tests may be indicated to differentiate
    - ▶ lung compliance
    - ▶ maximal inspiratory and expiratory pressure
  - ◆ neuromuscular conditions
    - ▶ normal compliance
    - ▶ decreased  $PI_{MAX}$ ,  $PE_{MAX}$

### Abnormal volumes - implications

- △ Reduced volumes - decreased compliance of lung and thorax ( $C_{LT}$ )
  - ◆ extrinsic restrictive disease
    - ▶ morbid obesity
    - ▶ pleural disease
    - ▶ thoracic deformity; e.g., kyphoscoliosis

### Abnormal volumes - implications

- △ Increased volumes - RV, FRC, TLC
  - ◆ emphysema - loss of lung elasticity
  - ◆ asthma (severe) - airway obstruction causing gas trapping

### Abnormal volumes - implications


- △ Reduced volumes
  - ◆ reduced effort - normal compliance ( $C_{LT}$ )
    - ▶ neuromuscular weakness; e.g., myasthenia gravis, muscular dystrophy may progress to extrinsic restrictive disease
    - ▶ purposeful - malingering (will not affect RV)

### Lung compliance

- △ Purpose - to determine restriction
- △ Measurement requires
  - ◆ static volumes (zero flow)
  - ◆ pleural pressure (esophageal pressure)
  - ◆ alveolar pressure (mouth pressure)
  - ◆ shutter for flow interruption
- △ Volume-pressure curve generated from data

## Lung compliance

### ▲ pleural pressure (esophageal pressure)



## Inspiratory/expiratory pressures

### ▲ $PI_{MAX}$ (NIF - NOT)

- ◆ maximal sustained pressure to generate inspiration
- ◆ measures strength of inspiratory muscles

### ▲ $PE_{MAX}$

- ◆ maximal sustained expiratory pressure
- ◆ reflects coughing capability

## Lung compliance

### ▲ Normals

- ◆  $C_L = 200 \text{ mL/cm H}_2\text{O}$
- ◆  $C_{TH} = 200 \text{ mL/cm H}_2\text{O}$
- ◆  $C_{LT} = 100 \text{ mL/cm H}_2\text{O}$

$$\begin{aligned} \text{▲ } 1/C_{LT} &= 1/C_L + 1/C_{TH} \implies \\ 1/100 &= 1/200 + 1/200 \end{aligned}$$

## Inspiratory/expiratory pressures

### ▲ Measurement devices

- ◆ pressure transducer OR
- ◆ aneroid manometer ( $\pm 200 \text{ cm H}_2\text{O}$ )
- ◆ mouthpiece (flanged) or mask
- ◆ nose clip

## Lung compliance

### ▲ Implications of abnormal compliance

- ◆ decreased  $C_L$  - intrinsic restrictive disease
- ◆ decreased  $C_{TH}$  - extrinsic restrictive disease
- ◆ increased  $C_L$  - emphysema

## Inspiratory/expiratory pressures

### ▲ Measurement technique

- ◆  $PI_{MAX}$ 
  - ▶ patient exhales to RV
  - ▶ airway occluded
  - ▶ maintain maximal pressure 2-3 s
- ◆  $PE_{MAX}$ 
  - ▶ patient inhales to TLC
  - ▶ airway occluded
  - ▶ maintain maximal pressure 2-3 s

## Inspiratory/expiratory pressures

### ^ Sources of error

- ◆ non-standardization of technique among operators
- ◆ inadequate cooperation
- ◆ leaks
- ◆ pressure generated by cheek muscles (mouthpiece measurement)

## Sniff nasal inspiratory pressure (SNIP)

### ^ Advantages

- ◆ natural maneuver - easy to do
- ◆ eliminates cheek muscle pressure
- ◆ reliability better than  $PI_{MAX}$
- ◆ validated for all age groups

## Normal values\*

Females	Mean (cm H <sub>2</sub> O)	SD (cm H <sub>2</sub> O)
$PI_{MAX}$	79	$\pm 19$
$PE_{MAX}$	111	$\pm 25$
Males		
$PI_{MAX}$	117	$\pm 25$
$PE_{MAX}$	192	$\pm 42$

\*Charfi MR, et al. 1991

## Pathophysiologic Basis For Spirometry

## Sniff nasal inspiratory pressure (SNIP)

### ^ Description

- ◆ manometer connected via tube to one nostril
- ◆ other nostril remains open
- ◆ patient instructed to sniff maximally

Assessing inspiratory stage of cough



## Physiology of air flow

### ^ Factors affecting the magnitude of airflow

- ◆ ventilatory muscle strength
- ◆ airway resistance
- ◆ elastic recoil of lung
- ◆ airway collapsibility

### Physiology of air flow

- ^ **Ventilatory muscle strength**
  - ◆ **inspiratory strength** - determines maximal inspired lung volume
    - ▶ greater volume ==> greater elastic recoil
    - ▶ greater volume ==> larger airways at beginning of expiration
  - ◆ **expiratory strength** - determines drive pressure for forced expiration

### Physiology of air flow

- ^ **Airway resistance ( $R_{AW}$ )** - pressure required to produce flow
  - ◆ predominant factor in  $R_{AW}$  - the radius of the airway
  - ◆  $R_{AW} \propto 1/r^4$  ==>
  - ◆ decreasing radius by 1/2 increases  $R_{AW}$  by  $2^4 = 16$  times ==>
  - ◆ pressure required to produce flow must increase by factor of 16 OR
  - ◆ the same pressure will produce 1/16 flow

### Physiology of air flow

- ^ **Airway resistance ( $R_{AW}$ )** - pressure required to produce flow
  - ◆ unit of measurement - cm H<sub>2</sub>O/L/sec
  - ◆ alternate unit - reciprocal of  $R_{AW}$  is airway conductance ( $G_{AW}$ )
  - ◆ measured directly with whole body plethysmography

### Physiology of air flow

- ^  $R_{AW}$  is non-constant during ventilatory cycles because the airway diameters change
  - ◆ high volumes - larger airways and lesser  $R_{AW}$
  - ◆ low volumes - smaller airways and greater  $R_{AW}$

### Physiology of air flow

- ^ **Airway resistance ( $R_{AW}$ )** - pressure required to produce flow
  - ◆ predominant factor in  $R_{AW}$  - the radius of the airway
  - ◆  $R_{AW} \propto 1/r^4$  ==>
  - ◆ decreasing radius by 1/2 increases  $R_{AW}$  by  $2^4 = 16$  times ==>

### Physiology of air flow

- ^ **Pathologic factors that affect  $R_{AW}$** 
  - ◆ bronchospasm - smooth muscle contraction
  - ◆ bronchiolitis - narrows lumens
  - ◆ airway edema
  - ◆ smooth muscle hypertrophy


## Physiology of air flow

### ▲ Pathologic factors that affect $R_{AW}$

- ◆ bronchospasm - smooth muscle contraction
- ◆ bronchiolitis - narrows lumens
- ◆ airway edema
- ◆ smooth muscle hypertrophy
- ◆ connective tissue remodeling
- ◆ increased numbers of mucus glands and goblet cells
- ◆ luminal secretions

## Physiology of air flow

### ▲ Migration of EPP



## Physiology of air flow


### ▲ Factors affecting expiratory airflow

- ◆ elastic recoil of lung - produces expiratory drive pressure during quiet breathing
- ◆ when elastic recoil is decreased and time for expiration is decreased, gas trapping occurs

## Physiology of air flow

### ▲ Inspiration vs. expiration


- ◆ during inspiration intrathoracic airways tend to expand and extrathoracic airways tend to collapse



## Physiology of air flow

### ▲ Airway collapsibility - migration of equal pressure point (EPP)

### ▲ EPP - point where intra-airway (drive) pressure equals surrounding pressure



## Physiology of air flow

### ▲ Fixed, vs. variable obstruction

- ◆ Fixed obstructions limit both inspiratory and expiratory flow; e.g.,:
  - ▶ tumors
  - ▶ congenital or acquired airway defects

## Physiology of air flow

- ^ Fixed, vs. variable obstruction
  - ◆ Fixed obstructions limit both inspiratory and expiratory flow; e.g.:
    - tumors
    - congenital or acquired airway defects
  - ◆ Variable obstructions change with intrathoracic pressure changes and vary with:
    - location of obstruction
    - airway collapsibility

## Spirometry indications

- ^ Patients presenting with undiagnosed respiratory symptoms; e.g.,
  - ◆ dyspnea
  - ◆ wheeze
  - ◆ cough

## Spirometry Methods and Devices

## Spirometry indications

- ^ Patients with suspected COPD, in particular those with a positive smoking history and:
  - ◆ increasing age
  - ◆ chronic cough
  - ◆ breathlessness on exertion, and
  - ◆ daily wheezing

## Spirometry

- ^ Description - measurement of inhaled and/or exhaled air flow
- ^ Purpose - to diagnose and evaluate airway obstructive disease
- ^ Types:
  - ◆ volume-time curves
  - ◆ flow-volume curves

## Spirometry indications

- ^ History of recurrent winter chest infections
- ^ Diagnosis of COPD
- ^ Monitoring patients with established obstructive disease
- ^ Diagnosis of asthma - normal spirometry results do NOT rule out asthma

### **Spirometry relative contraindications**

- ^ Known or suspected respiratory infection
- ^ Hemoptysis of unknown origin
- ^ Pneumothorax
- ^ Uncontrolled hypertension or history of hemorrhagic CVA

### **Spirometry devices**

- ^ Displays - screen, paper
- ^ Data recording - paper, digital
- ^ Data reporting - printout, digital file
- ^ Data storage - preferably digital

### **Spirometry relative contraindications**

- ^ Known or suspected respiratory infection
- ^ Hemoptysis of unknown origin
- ^ Pneumothorax
- ^ Uncontrolled hypertension or history of hemorrhagic CVA
- ^ Recent thoracic, abdominal or eye surgery
- ^ Nausea, vomiting, pain
- ^ Confusion, dementia

### **Spirometry procedures**

- ^ Pre-testing history and physical examination
- ^ Ask time of most recent bronchodilator
- ^ Measure height, weight

### **Spirometry devices**

- ^ Volume displacement - laboratories
- ^ Turbine spirometers
- ^ Pneumotachographs
- ^ Ultrasonic types
- ^ Heated wire type
- ^ Pitot tube type

### **Spirometry procedures**

- ^ Rapid inspiration, then forced expiration
- ^ Repeat three times
- ^ Two largest FVC, FEV<sub>1</sub> must be within 150 mL
- ^ Report greatest FVC, FEV<sub>1</sub>
- ^ Administer bronchodilator, if ordered and repeat FVC maneuvers


## Spirometry Interpretation

### FVC parameters

- △ **FVC\*** - forced vital capacity (liters)
- △ **FEV<sub>1</sub>\*** - forced exhaled volume in 1 s (liters)
- △ **FEV<sub>1</sub>/FVC%\***  $(FEV_1/FVC) * 100$

\*denotes undisputed utility

### Timed forced vital capacity (FVC)




Slope of curve at any point =  $dV/dT = \text{flow}$

AP Jones

### FVC parameters

- △ **FEV<sub>0.5</sub>** - forced exhaled volume in 0.5s (liters)
- △ **FEV<sub>6</sub>** - forced exhaled volume in 6s (liters)
- △ **FEV<sub>6</sub>/FVC%**

### Timed forced vital capacity (FVC)



AP Jones

### FVC parameters

- △ **FEV<sub>0.5</sub>** - forced exhaled volume in 0.5s (liters)
- △ **FEV<sub>6</sub>** - forced exhaled volume in 6s (liters)
- △ **FEV<sub>6</sub>/FVC%**
- △ **FEF<sub>25-75%</sub>** - forced expiratory flow between 25% and 75% expiration
- △ **PEF** - peak expiratory flow
- △ **V<sub>MAX50%</sub>, V<sub>MAX75%</sub>** - flows at 50% and 75% expiration respectively (flow-volume study)

## Interpretation

- ^ Assumed preconditions:
  - ◆ test is properly done by skilled operator
  - ◆ testing equipment is mechanically operational and calibrated
  - ◆ patient is capable, properly instructed and motivated

## Flow-Volume Studies

## Interpretation

- ^ FEV<sub>1</sub>/FVC is most reliable indicator for expiratory airway obstruction
- ^ FEF<sub>25-75%</sub> - its utility is controversial
- ^ FEV<sub>6</sub>/FVC - informing patients of lung age purported to aid smoking cessation??? (controversial)

## Flow-volume studies

- ^ Plot inspiratory and expiratory flows, versus lung volumes
- ^ Advantages:
  - ◆ display instantaneous flows at all points of FEVC and FIVC
  - ◆ graphic display of volume vs. flow for qualitative analysis
  - ◆ evaluates inspiratory and expiratory obstruction
  - ◆ capability to superimpose curves

## Interpretation

Degree of severity	FEV1% pred.
Mild	>80
Moderate	50-79
Severe	30-49
Very severe	<30


## Devices

- ^ Flow measuring device
- ^ Microprocessor, to:
  - ◆ process signals
  - ◆ generate display
  - ◆ store and recall data

## Technique


- ▲ Prepare patient
- ▲ Calibrate sensor
- ▲ Maximal inspiration, then,
- ▲ Forced expiration (FEVC), then,
- ▲ Forced inspiration (FIVC)

## Flow-volume loops




AP Jones

## Flow-volume loops




AP Jones

## Flow-volume loops



AP Jones

## Flow-volume loops



AP Jones

Variable obstruction will change with pressure  
Fixed obstruction remains constant

**Airways Resistance  
Studies**

### Airways resistance studies

- ▲ The gold standard for evaluating airway function
- ▲ Purposes:
  - ◆ asthma diagnosis
  - ◆ bronchodilator response
  - ◆ evaluate obstruction for patients who are unable or unwilling to generate FVC
  - ◆ bronchodilator response in patients who show a clinical response; but unchanged spirometry

### Parameters

- ▲  $R_{AW}$  - airways resistance (cm H<sub>2</sub>O/l/sec)
- ▲  $R_{AW0.5}$  - airways resistance at flow = 0.5 L/sec
- ▲  $sR_{AW}$  - specific airways resistance (cm H<sub>2</sub>O/L/sec/L)

### Plethysmographic measurement

- ▲ Advantages
  - ◆ effort independence
  - ◆ feasible with infants
  - ◆ rapid test completion
- ▲ Disadvantages
  - ◆ equipment expense
  - ◆ equipment space
  - ◆ extensive operator training

### Parameters

- ▲  $R_{AW}$  - airways resistance (cm H<sub>2</sub>O/l/sec)
- ▲  $R_{AW0.5}$  - airways resistance at flow = 0.5 L/sec
- ▲  $sR_{AW}$  - specific airways resistance (cm H<sub>2</sub>O/L/sec/L)
- ▲  $G_{AW}$  - conductance (L/sec/cm H<sub>2</sub>O)
- ▲  $sG_{AW}$  - specific conductance (L/sec/cm H<sub>2</sub>O/L)
- ▲ Shapes of loops

### Measurement

- ▲ Whole body plethysmographic measurement of:
  - ◆ alveolar pressure
  - ◆ gas flow
- ▲ Patient pants or breathes at TV
- ▲ Pressure-flow tracings generated

### Normal values

- ▲  $R_{AW0.5}$  - 0.6 - 2.4 cm H<sub>2</sub>O/L/sec
- ▲  $sG_{AW}$  - > 0.15 L/sec/cm H<sub>2</sub>O/L

## Summary & Review

- ^ Purposes (indications) for PFTs
- ^ Types of tests
- ^ Lung volume measurements
  - ◆ volumes & capacities
  - ◆ direct spirometry
  - ◆ direct spirometry devices
  - ◆ gas dilution methods for FRC
  - ◆ whole body plethysmography for TGV
  - ◆ alternative methods for volumes

## Summary & Review

- ^ Spirometry methods and devices
  - ◆ indications and contraindications for spirometry
  - ◆ spirometry devices
  - ◆ procedures
- ^ Spirometry interpretation
  - ◆ FVC curves
  - ◆ FVC parameters
  - ◆ Values for FEV<sub>1</sub>/FVC

## Summary & Review

- ^ Implications for abnormal volumes
- ^ Laboratory lung compliance measurement
- ^ Inspiratory and expiratory pressures

## Summary & Review

- ^ Flow-volume loops
  - ◆ restriction
  - ◆ variable intrathoracic obstruction
  - ◆ variable extrathoracic obstruction
  - ◆ fixed large airway obstruction
- ^ Airway resistance studies - gold standard
  - ◆ indications
  - ◆ procedure
  - ◆ parameters and interpretation

## Summary & Review

- ^ Pathophysiologic basis for spirometry - factors that affect airflow
  - ◆ muscle strength
  - ◆ airway resistance
  - ◆ elastic recoil of lung
  - ◆ airway collapsibility

## References

- ^ Ruppel GL. Manual of Pulmonary Function Testing, 9th Ed. Chaps. 2, 3, 5, 7, 9, 10. 2009 Mosby-Elsevier; St. Louis.
- ^ Goldman MD. Body plethysmography. The buyer's guide to respiratory care products.  
[http://dev.ersnet.org/uploads/Document/53/WEB\\_CHEMIN\\_2555\\_1194522257.pdf](http://dev.ersnet.org/uploads/Document/53/WEB_CHEMIN_2555_1194522257.pdf)
- ^ Levy ML, Quanjer PH, Booker R, Cooper BG, Holmes S, Small I. Diagnostic Spirometry in Primary Care: Proposed standards for general practice compliant with American Thoracic Society and European Respiratory Society recommendations. Prim Care Resp J 2009;18:130-147.
- ^ Pride, N. B. Contributions of Loss of Lung Recoil and of Enhanced Airways Collapsibility to the Airflow Obstruction of Chronic Bronchitis and Emphysema. J Clin Invest 1973; 52:2117-2128.

### References

- ^ Postma DS, Kerstjens HAM. Characteristics of Airway Hyperresponsiveness in Asthma and Chronic Obstructive Pulmonary Disease. *Am. J. Respir. Crit. Care Med.* 1998;158: S187-192S.
- ^ Cosio Piqueras, M.G., Cosio, M.G. Disease of the airways in chronic obstructive pulmonary disease. *Eur Respir J* 2001 18: 41S-49
- ^ MacIntyre NR, Selecky PA. Is there a role for screening spirometry? *Respir Care.* 2010 Jan;55(1):35-42.
- ^ Hansen JE, Sun XG, Wasserman K. Calculating gambling odds and lung ages for smokers. *Eur Respir J.* 2009 Sep 24.
- ^ Jing JY, Huang TC, Cui W, Xu F, Shen HH. Should FEV1/FEV6 replace FEV1/FVC ratio to detect airway obstruction? A metaanalysis. *Chest.* 2009 Apr;135(4):991-8

### References

- ^ Mayos M, Giner J, Casan P, Sanchis J. Measurement of maximal static respiratory pressures at the mouth with different air leaks. *Chest.* 1991 Aug;100(2):364-6.
- ^ Charfi MR, Matran R, Regnard J, Richard MO, Champeau J, Dall'ava J, Lockhart A. Maximal ventilatory pressure through the mouth in adults: normal values and explanatory variables. *Rev Mal Respir.* 1991;8(4):367-74.
- ^ Wanger JL, et al. Standardization of the measurement of lung volumes. *Eur Respir J* 2005;26:511-522.
- ^ Oppenheimer BW, Goldring RM, Herberg ME, Hofer IS, Reyfman PA, Liataud S, Rom WN, Reibman J, Berger KI. Distal airway function in symptomatic subjects with normal spirometry following World Trade Center dust exposure. *Chest.* 2007 Oct;132(4):1275-82. Epub 2007 Sep 21.

### References

- ^ Toda R, Hoshino T, Kawayama T, Imaoka H, Sakazaki Y, Tsuda T, Takada S, Kinoshita M, Iwanaga T, Aizawa H. Validation of "lung age" measured by spirometry and handy electronic FEV1/FEV6 meter in pulmonary diseases. *Intern Med.* 2009;48(7):513-21.
- ^ Schneider A, Gindner L, Tilemann L, Schermer T, Dinant GJ, Meyer FJ, Szecsenyi J. Diagnostic accuracy of spirometry in primary care. *BMC Pulm Med.* 2009 Jul 10;9:31.
- ^ O'Donnell CR, Bankier AA, Stiebellehner L, Reilly JJ, Brown R, Loring SH. Comparison of Plethysmographic and Helium Dilution Lung Volumes: Which is Best in COPD? *Chest.* 2009 Dec 18.
- ^ Hayes D Jr, Kraman SS. The physiologic basis of spirometry. *Respir Care.* 2009 Dec;54(12):1717-26.

### References

- ^ Brzozowska A, Majak P, Grzelewski T, Stelmach W, Kaczmarek J, Stelmach P, Jerzynska J, Stelmach I. Measurement of specific airway resistance decreased the risk of delay in asthma diagnosis in children. *Allergy Asthma Proc.* 2009 Jan-Feb;30(1):47-54.
- ^ Stocks J, Godfrey S, Beardsmore C, Bar-Yishay E, Castile R. Plethysmographic measurements of lung volume and airway resistance *Eur Respir J* 2001 17: 302-312.
- ^ Pakhale S, Bshouty Z, Marras TK. Comparison of per cent predicted and percentile values for pulmonary function test interpretation. *Can Respir J.* 2009 Nov-Dec;16(6):189-93.

### References

- ^ Brown MS, Kim HJ, Abtin F, Da Costa I, Pais R, Ahmad S, Angel E, Ni C, Kleerup EC, Gjertson DW, McNitt-Gray MF, Goldin JG. Reproducibility of Lung and Lobar Volume Measurements Using Computed Tomography. *Acad Radiol.*
- ^ Lee JS, Ra SW, Chae EJ, Seo JB, Lim SY, Kim TH, Lee SD, Oh YM. Validation of the Lower Limit of Normal Diffusing Capacity for Detecting Emphysema. *Respiration.* 2010 Jan 29.
- ^ Liptay MJ, Basu S, Hoaglin MC, Freedman N, Faber LP, Warren WH, Hammoud ZT, Kim AW. Diffusion lung capacity for carbon monoxide (DLCO) is an independent prognostic factor for long-term survival after curative lung resection for cancer. *J Surg Oncol.* 2009 Dec 15;100(8):703-7.